

## FINSEN SELFRETAINING RETRACTOR

blunt; 3X4 teeth

7.0cm (2.8")




**SKU:** 070-611-07

**Categories:** [Open Surgery](#), [Retractors](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>Aesculap</b>	BV003
<b>KLS Martin</b>	15-703-07
<b>AUTHORS</b>	<a href="#">FINSEN</a>
<b>Medicon</b>	24.04.08
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">toothed</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Selfretaining Retractors</a>
<b>Aesculap</b>	BV003
<b>KLS Martin</b>	15-703-07

## AUTHORS

[FINSEN](#)

## Medicon

24.04.08

## SPECIFICATIONS

[blunt](#), [toothed](#)