

FISCHER TENDON STRIPPER

Ø 4.0mm; straight

31.0cm (12.4")

SKU: 195-262-31

Categories: [Open Surgery](#), [Orthopaedics](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
KLS Martin	23-660-31
AUTHORS	FISCHER
Medicon	78.22.30
SPECIFICATIONS	straight

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Tendon Stripper
KLS Martin	23-660-31
AUTHORS	FISCHER
Medicon	78.22.30

SPECIFICATIONS

[straight](#)