

## FLEXOMETER

measuring range 150mm

23.5cm (9.4")

**SKU:** 100-235-23

**Categories:** [Open Surgery](#), [Diagnostics](#)



## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Measuring instruments</a>
<b>Aesculap</b>	AA846R
<b>Medicon</b>	44.38.05

## ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Measuring instruments</a>
<b>Aesculap</b>	AA846R
<b>Medicon</b>	44.38.05