

FLUTED BIT ELECTRODE; BIPOAR; INSERT

bit 1mm; for Ø 5.0mm shaft

for 628-201-53 and 628-205-53

SKU: 628-242-53

Categories: [Endoscopy](#), [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | HF-Instruments, bipolar |
| TYPE | Electrodes |

ZUSÄTZLICHE INFORMATIONEN

| | |
|--------------|---|
| UNIT | pcs. |
| RANGE | HF-Instruments, bipolar |
| TYPE | Electrodes |