

FLUTED BIT ELECTRODE BIPOLAR; 1MM; INSERT

short insulated; for coagulation forceps

628-101-00 and 628-105-00; 330mm

No CE mark! No labelling!



SKU: 628-142-53

Categories: [Endoscopy](#), [Laparoscopy](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Instruments, bipolar
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Instruments, bipolar
TYPE	Electrodes