

## FOERSTER-BALLENGER SPONGE FORCEPS

straight; serrated

25.0cm (10")


**SKU:** 030-160-25

**Categories:** [Open Surgery](#), [Hemostats](#)



### ZUSÄTZLICHE INFORMATIONEN

|                       |                                                      |
|-----------------------|------------------------------------------------------|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                 |
| <b>TYPE</b>           | <a href="#">Sponge and Dressing Forceps</a>          |
| <b>Aesculap</b>       | BF122                                                |
| <b>KLS Martin</b>     | 14-352-25                                            |
| <b>AUTHORS</b>        | <a href="#">BALLENGER</a> , <a href="#">FOERSTER</a> |
| <b>Medicon</b>        | 19.50.24                                             |
| <b>SPECIFICATIONS</b> | <a href="#">serrated</a> , <a href="#">straight</a>  |



### ZUSÄTZLICHE INFORMATIONEN

|                   |                                             |
|-------------------|---------------------------------------------|
| <b>UNIT</b>       | <a href="#">pcs.</a>                        |
| <b>TYPE</b>       | <a href="#">Sponge and Dressing Forceps</a> |
| <b>Aesculap</b>   | BF122                                       |
| <b>KLS Martin</b> | 14-352-25                                   |

## AUTHORS

[BALLENGER](#), [FOERSTER](#)

## Medicon

19.50.24

## SPECIFICATIONS

[serrated](#), [straight](#)