

**FOERSTER-BALLENGER SPONGE FORCEPS**

straight; serrated


18.0cm (7.2")



**SKU:** 030-160-18  
**Categories:** [Open Surgery](#), [Hemostats](#)

**ZUSÄTZLICHE INFORMATIONEN**

|                |                                                      |
|----------------|------------------------------------------------------|
| UNIT           | <a href="#">pcs.</a>                                 |
| TYPE           | <a href="#">Sponge and Dressing Forceps</a>          |
| Aesculap       | BF112                                                |
| KLS Martin     | 14-352-18                                            |
| AUTHORS        | <a href="#">BALLENGER</a> , <a href="#">FOERSTER</a> |
| Medicon        | 19.50.18                                             |
| SPECIFICATIONS | <a href="#">serrated</a> , <a href="#">straight</a>  |



**ZUSÄTZLICHE INFORMATIONEN**

|            |                                             |
|------------|---------------------------------------------|
| UNIT       | <a href="#">pcs.</a>                        |
| TYPE       | <a href="#">Sponge and Dressing Forceps</a> |
| Aesculap   | BF112                                       |
| KLS Martin | 14-352-18                                   |

## AUTHORS

[BALLENGER](#), [FOERSTER](#)

## Medicon

19.50.18

## SPECIFICATIONS

[serrated](#), [straight](#)