

FOMON ALAR HOOKLET

2-pronged; button ended

16.0cm (6.4")


SKU: 112-305-16

Categories: [Open Surgery](#), [Rhinology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Delicate Hooks
Aesculap	OL202
KLS Martin	37-706-16
AUTHORS	FOMON
SPECIFICATIONS	toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Delicate Hooks
Aesculap	OL202
KLS Martin	37-706-16
AUTHORS	FOMON

SPECIFICATIONS

[toothed](#)