

## FOMON ALAR HOOKLET

button ended; blunt

17.0cm (6.8")

**SKU:** 112-303-17

**Categories:** [Open Surgery](#), [Rhinology](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>Aesculap</b>	VF225
<b>KLS Martin</b>	37-694-17
<b>AUTHORS</b>	<a href="#">FOMON</a>
<b>SPECIFICATIONS</b>	<a href="#">blunt</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Delicate Hooks</a>
<b>Aesculap</b>	VF225
<b>KLS Martin</b>	37-694-17
<b>AUTHORS</b>	<a href="#">FOMON</a>

## SPECIFICATIONS

[blunt](#)