

## FOMON CHISEL

straight; 7.0mm

16.0cm (6.4")

**SKU:** 112-537-07

**Categories:** [Open Surgery](#), [Rhinology](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Chisels and Gouges</a>
<b>Aesculap</b>	OL339
<b>KLS Martin</b>	37-623-07
<b>AUTHORS</b>	<a href="#">FOMON</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Chisels and Gouges</a>
<b>Aesculap</b>	OL339
<b>KLS Martin</b>	37-623-07
<b>AUTHORS</b>	<a href="#">FOMON</a>