

FOREIGN BODY FORCEPS; Ø 1.2 MM; FL 600MM

4 Fr.; spoon serrated fine;long jaw part

semi-rigid;double action; with LuerLock

insert only



SKU: 583-472-12LLA
Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscopy
TYPE	Semi-rigid Instruments



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscopy
TYPE	Semi-rigid Instruments