

FOREIGN BODY FORCEPS; Ø 1.7MM; FL 600MM

spoon serrated; long jaw part;semi-rigid

double action; with LuerLock;with handle



with Luer Lock

SKU: 583-472-16LL

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscapy
TYPE	Semi-rigid Instruments



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	Uretero-Renoscapy
TYPE	Semi-rigid Instruments