

## FOREIGN BODY FORCEPS; Ø 1.7MM; FL 600MM

5 Fr.; spoon serrated; long jaw part

semi-rigid; double action; with handle



**SKU:** 583-472-16

**Categories:** [Endoscopy](#), [Urology](#)

### ZUSÄTZLICHE INFORMATIONEN

|              |                                        |
|--------------|----------------------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>                   |
| <b>RANGE</b> | <a href="#">Uretero-Renoscapy</a>      |
| <b>TYPE</b>  | <a href="#">Semi-rigid Instruments</a> |



### ZUSÄTZLICHE INFORMATIONEN

|              |                                        |
|--------------|----------------------------------------|
| <b>UNIT</b>  | <a href="#">pcs.</a>                   |
| <b>RANGE</b> | <a href="#">Uretero-Renoscapy</a>      |
| <b>TYPE</b>  | <a href="#">Semi-rigid Instruments</a> |