

FORIGN BODY FORCEPS; Ø 1.7MM; FL 600MM

spoon serrated; long jaw part;semi-rigid

double action; with LuerLock;insert only



with Luer Lock

SKU: 583-472-16LLA

Categories: [Endoscopy](#), [Urology](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT

[pcs.](#)

RANGE

[Uretero-Renoscapy](#)

TYPE

[Semi-rigid Instruments](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT

[pcs.](#)

RANGE

[Uretero-Renoscapy](#)

TYPE

[Semi-rigid Instruments](#)