

FRAZIER SUCTION TUBE

15 Fr.; 45° angled

10.0cm (4") working length

SKU: 085-030-15

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| Aesculap | GF925 |
| KLS Martin | 18-523-13 |
| AUTHORS | FRAZIER |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| Aesculap | GF925 |
| KLS Martin | 18-523-13 |
| AUTHORS | FRAZIER |