

## FRAZIER SUCTION TUBE

15 Fr.; 45° angled

13.0cm (5.2") working length

**SKU:** 085-031-15

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)



### ZUSÄTZLICHE INFORMATIONEN

|                |                               |
|----------------|-------------------------------|
| <b>UNIT</b>    | <a href="#">pcs.</a>          |
| <b>TYPE</b>    | <a href="#">Suction Tubes</a> |
| <b>AUTHORS</b> | <a href="#">FRAZIER</a>       |
| <b>Medicon</b> | 09.83.34                      |

### ZUSÄTZLICHE INFORMATIONEN

|                |                               |
|----------------|-------------------------------|
| <b>UNIT</b>    | <a href="#">pcs.</a>          |
| <b>TYPE</b>    | <a href="#">Suction Tubes</a> |
| <b>AUTHORS</b> | <a href="#">FRAZIER</a>       |
| <b>Medicon</b> | 09.83.34                      |