

## FRAZIER SUCTION TUBE

3 Fr.; 45° angled

18.0cm (7.2") working length



**SKU:** 085-033-03

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)

### ZUSÄTZLICHE INFORMATIONEN

|                |                               |
|----------------|-------------------------------|
| <b>UNIT</b>    | <a href="#">pcs.</a>          |
| <b>TYPE</b>    | <a href="#">Suction Tubes</a> |
| <b>AUTHORS</b> | <a href="#">FRAZIER</a>       |

### ZUSÄTZLICHE INFORMATIONEN

|                |                               |
|----------------|-------------------------------|
| <b>UNIT</b>    | <a href="#">pcs.</a>          |
| <b>TYPE</b>    | <a href="#">Suction Tubes</a> |
| <b>AUTHORS</b> | <a href="#">FRAZIER</a>       |