

FRAZIER SUCTION TUBE

4 Fr.; 45° angled

18.0cm (7.2") working length

SKU: 085-033-04

Categories: [Open Surgery](#), [Puncture / Special Needles](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|----------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| AUTHORS | FRAZIER |

ZUSÄTZLICHE INFORMATIONEN

| | |
|----------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Suction Tubes |
| AUTHORS | FRAZIER |