

FRAZIER SUCTION TUBE

6 Fr.; 45° angled

18.0cm (7.2") working length



SKU: 085-033-06

Categories: [Open Surgery](#), [Puncture / Special Needles](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
KLS Martin	24-763-17
AUTHORS	FRAZIER
Medicon	09.83.76

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
KLS Martin	24-763-17
AUTHORS	FRAZIER
Medicon	09.83.76