

## FRAZIER SUCTION TUBE

8 Fr.; 45° angled

10.0cm (4") working length



**SKU:** 085-030-08

**Categories:** [Open Surgery](#), [Puncture / Special Needles](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>Aesculap</b>	GF928
<b>KLS Martin</b>	18-523-08
<b>AUTHORS</b>	<a href="#">FRAZIER</a>
<b>Medicon</b>	09.83.48

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Suction Tubes</a>
<b>Aesculap</b>	GF928
<b>KLS Martin</b>	18-523-08
<b>AUTHORS</b>	<a href="#">FRAZIER</a>

**Medicon**

09.83.48