

GIL-VERNET RETRACTOR

23mm; flexible


24.0cm (9.6")

SKU: 062-230-23
Categories: [Open Surgery](#), [Retractors](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Delicate Hooks
KLS Martin	15-174-23
AUTHORS	GIL-VERNET



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Delicate Hooks
KLS Martin	15-174-23
AUTHORS	GIL-VERNET