

GLOVER BULLDOG CLAMP

straight; atraumatic

9.5cm (3.7")

SKU: 026-298-09

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB434
KLS Martin	24-398-09
AUTHORS	GLOVER
Medicon	15.18.09
SPECIFICATIONS	atraumatic , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	FB434
KLS Martin	24-398-09

AUTHORS

[GLOVER](#)

Medicon

15.18.09

SPECIFICATIONS

[atraumatic](#), [straight](#)