

GRASPING FORCEP ELECTRODE;BIPOLAR;INSERT

fenestrated; for Ø 5.0mm shaft
 for 628-201-53 and 628-205-53

SKU: 628-244-53
Categories: [Endoscopy](#), [Laparoscopy](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Instruments, bipolar
TYPE	Electrodes

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
RANGE	HF-Instruments, bipolar
TYPE	Electrodes