

GUILLEN SUCTION RASPATORY

straight tip; Ø 4.0mm; 20.0cm



SKU: 085-074-20
Categories: [Open Surgery](#), [Puncture / Special Needles](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
KLS Martin	37-840-20
AUTHORS	GUILLEN

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Suction Tubes
KLS Martin	37-840-20
AUTHORS	GUILLEN