

GUYON KIDNEY PEDICLE CLAMP

curved; serrated

24.0cm (9.6")



SKU: 040-509-24

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Kidney Pedicle Clamp Forceps
Aesculap	BH839
KLS Martin	13-531-24
AUTHORS	GUYON
Medicon	51.56.24
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Kidney Pedicle Clamp Forceps
Aesculap	BH839
KLS Martin	13-531-24

AUTHORS

[GUYON](#)

Medicon

51.56.24

SPECIFICATIONS

[curved](#), [serrated](#)