

## GUYON KIDNEY PEDICLE CLAMP

curved; serrated

23.0cm (9.2")


**SKU:** 040-513-23

**Categories:** [Open Surgery](#), [Hemostats](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Kidney Pedicle Clamp Forceps</a>
<b>AUTHORS</b>	<a href="#">GUYON</a>
<b>Medicon</b>	51.56.15
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Kidney Pedicle Clamp Forceps</a>
<b>AUTHORS</b>	<a href="#">GUYON</a>
<b>Medicon</b>	51.56.15
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>