

GUYON KIDNEY PEDICLE CLAMP

curved; serrated

20.0cm (8")



SKU: 040-509-20

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Kidney Pedicle Clamp Forceps
Aesculap	BH835
KLS Martin	13-531-20
AUTHORS	GUYON
Medicon	51.56.20
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Kidney Pedicle Clamp Forceps
Aesculap	BH835
KLS Martin	13-531-20

AUTHORS

[GUYON](#)

Medicon

51.56.20

SPECIFICATIONS

[curved](#), [serrated](#)