

GUYON PEDICLE CLAMP

atraumatic; curved

23.0cm (9.2")




SKU: 040-519-23

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Kidney Pedicle Clamp Forceps |
| Aesculap | EF002 |
| KLS Martin | 13-529-24 |
| AUTHORS | GUYON |
| Medicon | 51.56.26 |
| SPECIFICATIONS | curved |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Kidney Pedicle Clamp Forceps |
| Aesculap | EF002 |
| KLS Martin | 13-529-24 |

AUTHORS

[GUYON](#)

Medicon

51.56.26

SPECIFICATIONS

[curved](#)