

HALSEY NEEDLE HOLDER

straight; serrated

13.0cm (5.2")



SKU: 050-040-13

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM205
AUTHORS	HALSEY
Medicon	10.12.13
SPECIFICATIONS	serrated , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM205
AUTHORS	HALSEY
Medicon	10.12.13

SPECIFICATIONS

[serrated](#), [straight](#)