

HALSTED HAEMOSTATIC FORCEPS

curved; serrated

21.0cm (8.4")




SKU: 040-051-21

Categories: [Open Surgery](#), [Hemostats](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Artery Forceps |
| Aesculap | BH211 |
| KLS Martin | 13-317-21 |
| AUTHORS | HALSTED |
| Medicon | 16.41.21 |
| SPECIFICATIONS | curved , serrated |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--------------------------------|
| UNIT | pcs. |
| TYPE | Artery Forceps |
| Aesculap | BH211 |
| KLS Martin | 13-317-21 |

AUTHORS

[HALSTED](#)

Medicon

16.41.21

SPECIFICATIONS

[curved](#), [serrated](#)