

HD MINI COMPACT NEPHROSCOPE WITH

lateral view, view 6°, 40.000 pixel

working length 220mm

working channel Ø 2.2mm, shaft Ø 4.0mm


SKU: 601-145-00

Categories: [Endoscopy](#), [Equipment](#), [Urology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Telescopes](#)
SCOPE [Urology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)
RANGE [Telescopes](#)
SCOPE [Urology](#)