

## HD MINI COMPACT NEPHROSCOPE WITH

lateral view, view 6°, 40.000 pixel

working length 220mm

working channel Ø 2.2mm, shaft Ø 4.0mm


**SKU:** 601-145-00

**Categories:** [Endoscopy](#), [Equipment](#), [Urology](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Telescopes</a>
<b>SCOPE</b>	<a href="#">Urology</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>RANGE</b>	<a href="#">Telescopes</a>
<b>SCOPE</b>	<a href="#">Urology</a>