

HEANEY NEEDLE HOLDER

curved; serrated

21.0cm (8.4")



SKU: 050-141-21

Categories: [Open Surgery](#), [Suture](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM259
KLS Martin	20-403-21
AUTHORS	HEANEY
Medicon	10.31.21
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Needle Holders
Aesculap	BM259
KLS Martin	20-403-21

AUTHORS

[HEANEY](#)

Medicon

10.31.21

SPECIFICATIONS

[curved](#), [serrated](#)