

HEGAR UTERINE DILATOR

with sloped handle; Ø 14.0mm
 18.5cm (7.25")

SKU: 180-350-14
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM114
KLS Martin	32-400-14
AUTHORS	HEGAR
Medicon	52.30.14
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM114
KLS Martin	32-400-14

AUTHORS

[HEGAR](#)

Medicon

52.30.14

SPECIFICATIONS

[sloped handle](#)