

HEGAR UTERINE DILATOR

with sloped handle; Ø 2.0mm
 18.5cm (7.25")

SKU: 180-350-02
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM102
KLS Martin	32-400-02
AUTHORS	HEGAR
Medicon	52.30.02
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM102
KLS Martin	32-400-02

AUTHORS

[HEGAR](#)

Medicon

52.30.02

SPECIFICATIONS

[sloped handle](#)