

HEGAR UTERINE DILATOR

with sloped handle; Ø 22.0mm

18.5cm (7.25")

SKU: 180-350-22

Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM122
KLS Martin	32-400-22
AUTHORS	HEGAR
Medicon	52.30.22
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM122
KLS Martin	32-400-22

AUTHORS

[HEGAR](#)

Medicon

52.30.22

SPECIFICATIONS

[sloped handle](#)