

HEGAR UTERINE DILATOR

with sloped handle; Ø 25.0mm
 18.5cm (7.25")

SKU: 180-350-25
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM125
KLS Martin	32-400-25
AUTHORS	HEGAR
Medicon	52.30.25
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM125
KLS Martin	32-400-25

AUTHORS

[HEGAR](#)

Medicon

52.30.25

SPECIFICATIONS

[sloped handle](#)