

HEGAR UTERINE DILATOR

with sloped handle; Ø 26.0mm
 18.5cm (7.25")

SKU: 180-350-26
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM126
KLS Martin	32-400-26
AUTHORS	HEGAR
Medicon	52.30.26
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM126
KLS Martin	32-400-26

AUTHORS

[HEGAR](#)

Medicon

52.30.26

SPECIFICATIONS

[sloped handle](#)