

HEGAR UTERINE DILATOR

with sloped handle; Ø 27.0mm  
 18.5cm (7.25")

SKU: 180-350-27  
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

|                |                                  |
|----------------|----------------------------------|
| UNIT           | <a href="#">pcs.</a>             |
| TYPE           | <a href="#">Uterine Dilators</a> |
| KLS Martin     | 32-400-27                        |
| AUTHORS        | <a href="#">HEGAR</a>            |
| SPECIFICATIONS | <a href="#">sloped handle</a>    |

ZUSÄTZLICHE INFORMATIONEN

|                |                                  |
|----------------|----------------------------------|
| UNIT           | <a href="#">pcs.</a>             |
| TYPE           | <a href="#">Uterine Dilators</a> |
| KLS Martin     | 32-400-27                        |
| AUTHORS        | <a href="#">HEGAR</a>            |
| SPECIFICATIONS | <a href="#">sloped handle</a>    |