

HEGAR UTERINE DILATOR

with sloped handle; Ø 30.0mm
 18.5cm (7.25")

SKU: 180-350-30
Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-400-30
AUTHORS	HEGAR
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-400-30
AUTHORS	HEGAR
SPECIFICATIONS	sloped handle