

HEGAR UTERINE DILATOR

with sloped handle; Ø 7.5mm
 18.5cm (7.25")

SKU: 180-351-07
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM137
KLS Martin	32-401-07
AUTHORS	HEGAR
Medicon	52.31.07
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM137
KLS Martin	32-401-07

AUTHORS

[HEGAR](#)

Medicon

52.31.07

SPECIFICATIONS

[sloped handle](#)