

HEGAR UTERINE DILATOR

with sloped handle; Ø 8.5mm
 18.5cm (7.25")

SKU: 180-351-08
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM138
KLS Martin	32-401-08
AUTHORS	HEGAR
Medicon	52.31.08
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM138
KLS Martin	32-401-08

AUTHORS

[HEGAR](#)

Medicon

52.31.08

SPECIFICATIONS

[sloped handle](#)