

HEGAR UTERINE DILATOR

with sloped handle; Ø 9.5mm
 18.5cm (7.25")

SKU: 180-351-09
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM139
KLS Martin	32-401-09
AUTHORS	HEGAR
Medicon	52.31.09
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM139
KLS Martin	32-401-09

AUTHORS

[HEGAR](#)

Medicon

52.31.09

SPECIFICATIONS

[sloped handle](#)