

HEGAR UTERINE DILATOR

with sloped handle; Ø 4.0mm
 18.5cm (7.25")

SKU: 180-350-04
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM104
KLS Martin	32-400-04
AUTHORS	HEGAR
Medicon	52.30.04
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM104
KLS Martin	32-400-04

AUTHORS

[HEGAR](#)

Medicon

52.30.04

SPECIFICATIONS

[sloped handle](#)