

HEGAR UTERINE DILATOR

with sloped handle; Ø 17.5mm
 18.5cm (7.25")

SKU: 180-351-17
Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM147
KLS Martin	32-401-17
AUTHORS	HEGAR
Medicon	52.31.17
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM147
KLS Martin	32-401-17

AUTHORS

[HEGAR](#)

Medicon

52.31.17

SPECIFICATIONS

[sloped handle](#)