

HEGAR UTERINE DILATOR

with sloped handle; Ø 18.5mm
 18.5cm (7.25")

SKU: 180-351-18
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-401-18
AUTHORS	HEGAR
Medicon	52.31.18
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
KLS Martin	32-401-18
AUTHORS	HEGAR
Medicon	52.31.18

SPECIFICATIONS

[sloped handle](#)