

HEGAR UTERINE DILATOR

with sloped handle; Ø 5.0mm
 18.5cm (7.25")

SKU: 180-350-05
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM105
KLS Martin	32-400-05
AUTHORS	HEGAR
Medicon	52.30.05
SPECIFICATIONS	sloped handle

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM105
KLS Martin	32-400-05

AUTHORS	HEGAR
Medicon	52.30.05
SPECIFICATIONS	sloped handle