

## HEGAR UTERINE DILATOR

double ended; Ø 5.0 / 6.0mm

20.0cm (8")

**SKU:** 180-354-05

**Categories:** [Open Surgery](#), [Gynecology](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Uterine Dilators</a>
<b>Aesculap</b>	EM233
<b>KLS Martin</b>	32-430-06
<b>AUTHORS</b>	<a href="#">HEGAR</a>
<b>Medicon</b>	52.35.05
<b>SPECIFICATIONS</b>	<a href="#">double ended</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Uterine Dilators</a>
<b>Aesculap</b>	EM233
<b>KLS Martin</b>	32-430-06

## AUTHORS

[HEGAR](#)

## Medicon

52.35.05

## SPECIFICATIONS

[double ended](#)