

HEGAR UTERINE DILATOR

double ended; Ø 7.0 / 8.0mm

20.0cm (8")

SKU: 180-354-07
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM234
KLS Martin	32-430-08
AUTHORS	HEGAR
Medicon	52.35.07
SPECIFICATIONS	double ended

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM234
KLS Martin	32-430-08

AUTHORS

[HEGAR](#)

Medicon

52.35.07

SPECIFICATIONS

[double ended](#)