

HEGAR UTERINE DILATOR

double ended; Ø 13.0 / 14.0mm
 20.0cm (8")

SKU: 180-354-13
 Categories: [Open Surgery](#), [Gynecology](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM237
KLS Martin	32-430-14
AUTHORS	HEGAR
Medicon	52.35.13
SPECIFICATIONS	double ended

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Uterine Dilators
Aesculap	EM237
KLS Martin	32-430-14

AUTHORS

[HEGAR](#)

Medicon

52.35.13

SPECIFICATIONS

[double ended](#)